

Introduction

NGS basic concepts and short-reads QC



The diagram illustrates the process of variant calling from a biological sample to genotype inference. It is divided into two main sections: a detailed view of a diploid genetic variant and a flowchart of the sequencing process.

A diploid genetic variant

This section shows a blue gear icon representing a biological sample, which leads to a DNA double helix structure labeled "DNA extraction". The extracted DNA is then shown as a pair of homologous chromosomes, each with a different colored band representing a variant. The variant is labeled "DNA variant" and is shown in three states: "AA hom" (homozygous), "AG het" (heterozygous), and "GG hom" (homozygous).

Sequencing Process Flowchart


The flowchart shows the process from a "Biological sample" (represented by blue gears) to a "Pool of molecules" (represented by DNA double helices). This pool is then subjected to "random sampling" (represented by a grid of colored dots) and "NGS sequencing" (represented by a blue machine with a DNA helix on its screen). The resulting data is then processed for "Variant calling" (represented by a grid of colored bars) and "Inference of genotype" (represented by a computer monitor icon).

VAF (Variant Allele Frequency)

The VAF is shown as a bar chart with two bars: a green bar for "A (0.6)" and an orange bar for "G (0.4)". The x-axis is labeled "SNP Allelic Ratio (AR)" and ranges from 0.000 to 0.995. The y-axis is labeled "Number of SNPs" and ranges from 0 to 6000. The chart shows a distribution of SNPs with a peak at 0.5, labeled "Expected Ref/Alt AR", and two shaded regions labeled "Unexpected AR".

Variant calling

The variant calling results are shown as a grid of colored bars representing the sequence of bases (A, C, G, T) for each variant. The bases are labeled "A C G T" at the bottom.



NGS basic concepts

NGS workflow

Essential steps in NGS data generation
for short and long-reads

Next-Generation Sequencing

Massive parallel sequencing of DNA fragments

2nd generation

Short fragments (50-300bp)

Needs DNA fragmentation




3rd generation

Long fragments (10-100kb or more)

DNA molecules are directly used for sequencing



NGS short reads - library preparation



NGS short reads - library preparation

LONG READS

Template DNA



Fragmentation

Adapters

Clonal
amplification


Sequencing

PacBio



SHORT READS

Single index



Dual index
(unique or combinatorial)



xGen UDI-UMI adapter



- Flow cell binding sequence: Platform-specific sequences for library binding to instrument
- Sequencing primer sites: Binding sites for general sequencing primers
- Sample indexes: Short sequences specific to a given sample library
- Molecular index/barcode: Short sequence used to uniquely tag each molecule in a given sample library
- Insert: Target DNA or RNA fragment from a given sample library

NGS short reads - library preparation

LONG READS

Template DNA

Fragmentation



Adapters

Clonal
amplification

Sequencing

**NO ACTION
REQUIRED**

SHORT READS



iontorrent
by Thermo Fisher Scientific

NGS short reads - library preparation

Template DNA


Fragmentation

Adapters

Clonal
amplification

Sequencing

LONG READS




ZMW wells
Sites where sequencing takes place

Labelled nucleotides
All four dNTPs are labelled and available for incorporation

Modified polymerase
As a nucleotide is incorporated by the polymerase, a camera records the emitted light


PacBio output
A camera records the changing colours from all ZMWs; each colour change corresponds to one base

PacBio




SHORT READS


illumina




Nucleotide addition
Fluorophore-labelled, terminally blocked nucleotides hybridize to complementary base. Each cluster on a slide can incorporate a different base.



Imaging
Slides are imaged with either two or four laser channels. Each cluster emits a colour corresponding to the base incorporated during this cycle.



Cleavage
Fluorophores are cleaved and washed from flow cells and the 3'-OH group is regenerated. A new cycle begins with the addition of new nucleotides.



Semiconductor sequencing
As a base is incorporated, a single H⁺ ion is released, which is detected by a CMOS-ISFET sensor

ion torrent
by Thermo Fisher Scientific

NGS reads structure

Understand the structure of sequencing libraries and the resulting data

- adapters
- read configuration
- UMIs

NGS data - understand library structure

Single index



Dual index
(unique or combinatorial)



xGen UDI-UMI adapter



- Flow cell binding sequence: Platform-specific sequences for library binding to instrument
- Sequencing primer sites: Binding sites for general sequencing primers
- Sample indexes: Short sequences specific to a given sample library
- Molecular index/barcode: Short sequence used to uniquely tag each molecule in a given sample library
- Insert: Target DNA or RNA fragment from a given sample library

Library barcode(s) for multiplexing


Example of sequencing reads

Standard library preparation


Read contains the PCR primer

Example of sequencing reads

Amplicon targeted library





NGS data - single end and paired-end short reads




In **single end protocol**, each DNA molecule is sequenced once starting from a specific end.

In **paired-end protocol**, it is sequenced twice starting from opposing ends

Illumina paired-end sequencing



NGS data - understand UMIs



NGS data - library structure impacts downstream analysis

Reference genome

... AGTAGCTACGGATTTCGGAATCATCGCAATTCCTTAGCTTGAGGCACA ...

AGGTTTTCGGAATCATCGCAATTCCTTAGC

Aligned read

Residual adapter sequence can
create a false SNV



NGS data - library structure impacts downstream analysis



- A sequencing read derives from multiple molecular manipulations and contains specific elements besides the template DNA
- Knowing the elements and structure of sequencing reads is crucial for proper downstream processing
- Fail to remove adapters parts can results in increased error rate in variant calling
- When a library is generated by targeted amplification, no variant can be detected in the target primers sequences.

Data formats


Main file formats

- FASTQ: sequences

Sequences - FASTQ files

Sequence and base quality of each read

```
@A02059:47:HTV2VDSX7:1:1101:2284:1000 1:N:0:TGTAATCGAC+NGCGGTGATC
GCCTGACCTCTCCATCAACTGTTACCGAATGGTTCATCCATGTTGGGTTTGCCTAAATCACTTTACATCATTAGAGTTTGAA
+
:FFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFF
@A02059:47:HTV2VDSX7:1:1101:4743:1000 1:N:0:TGTAATCGAC+NGCGGTGATC
GGCCTCGGCCAGGCACGGTGGCTCATGCCGTGAATCCAGCACTTTGGGAGGCCGAGGCAGGCAGATCACCTGAGGTCAGGAGT
+
FFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFF
@A02059:47:HTV2VDSX7:1:1101:5195:1000 1:N:0:TGTAATCGAC+NGCGGTGATC
GTTGAAGTCTGGAACCTTTTAGGATCCTTTAAATGACTAAATGTTAATGTTTGGATAATTATAATGCTTAACATCTA
+
FFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFF
@A02059:47:HTV2VDSX7:1:1101:8323:1000 1:N:0:TGTAATCGAC+NGCGGTGATC
ACCTGACCCCTTGATGCCCTTATTCCTCTTCGCCTCATCTCCTCAGTGGGGTTATTAGTCTGATTCACTCTATCATTTCTCTAA
+
FFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFFF
```



- **Label**
unique identifier for the single read (instrument ID, run number, chip ID, tile, tile XY ...)
- **Q scores**
phred like quality score for each base sequenced encoded as ASCII character

Sequences - Base quality in FASTQ files

PHRED quality

$$Q = -10 \log_{10} P$$

$$P = 10^{-\frac{Q}{10}}$$

Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10,000	99.99%

Quality encoding as ASCII characters

```
@SRR038845.3 HWI-EAS038:6:1:0:1938 length=36  
CAACGAGTTACACCTTGGCCGACAGGCCCGGGTAA  
+SRR038845.3 HWI-EAS038:6:1:0:1938 length=36  
BA@7>B=>:>>7@7@>>9=BAA?;>52;>:9=8.=A
```

Quality converted to single ASCII character
PHRED+33 \Rightarrow ASCII code \Rightarrow Char

Char	ASCII code	Phred quality
B	66	33
A	65	32
@	64	31
7	55	22

ASCII Table

Dec	Hx	Oct	Html	Chr	Dec	Hx	Oct	Html	Chr
64	40	100	@	@	96	60	140	`	`
65	41	101	A	A	97	61	141	a	a
66	42	102	B	B	98	62	142	b	b
67	43	103	C	C	99	63	143	c	c
68	44	104	D	D	100	64	144	d	d
69	45	105	E	E	101	65	145	e	e
70	46	106	F	F	102	66	146	f	f
71	47	107	G	G	103	67	147	g	g
72	48	110	H	H	104	68	150	h	h
73	49	111	I	I	105	69	151	i	i
74	4A	112	J	J	106	6A	152	j	j
75	4B	113	K	K	107	6B	153	k	k
76	4C	114	L	L	108	6C	154	l	l
77	4D	115	M	M	109	6D	155	m	m
78	4E	116	N	N	110	6E	156	n	n
79	4F	117	O	O	111	6F	157	o	o
80	50	120	P	P	112	70	160	p	p
81	51	121	Q	Q	113	71	161	q	q
82	52	122	R	R	114	72	162	r	r

Short-reads QC

Reads cleaning


Remove unwanted sequences and poor quality bases

- adapter trimming
- quality trimming
- fixed length trimming

Useful tools

fastp, cutadapt

Adapter trimming - remove unwanted fixed sequences



Example from [cutadapt tool](#)

Options -a and -b

Read runs into adapter



Adapter

Removed sequence

Adapter within read



Option -a

Full adapter in the beginning



Option -b

Full adapter in the beginning




Partial adapter in the beginning






Trimming tools scan the read ends removing fixed sequences


Based on adapter known sequences the algorithm search for partial matches and then extend the matching region to completely remove adapter sequences



Quality trimming - remove low quality bases from read ends



Fixed length trimming - remove a fixed amount of bases from read ends



Short-reads QC

- raw reads length and quality
- GC content
- adapter content
- mapping statistics

Useful tools



fastp, fastQC, samtools

Reads QC - Per base sequence quality across the reads

- Ideally, base quality should be ≥ 30 across all the read
- A little decrease in quality is expected toward the end of the read
- Base quality is considered during variant calling and can affect variant caller performances





Per base sequence quality



Read QC - sequence length distribution

- Most read should have the expected read length
- A small tail on the left is acceptable, especially after trimming




*Sequencing problems,
Low quality reads
that got trimmed,
adapter dimers*



Reads QC - Per base sequence content across the reads

In most cases random sequences are expected
→ balanced base composition across the reads

Per base sequence content



Residual fixed elements or sequencing artefacts




Reads QC - reads GC content


- Distribution should follow the expected for the sequenced organism
- Peak around 45% for human genome samples



Per sequence GC content





Contamination from DNA of other organisms or sequencing artefacts




Read QC - residual adapter content

- Low presence of adapter sequences at the read ends
- Residual adapters can be cleaned by trimming



Aligned reads QC - mapping statistics



All libraries

- Mapping quality distribution
(≥ 30 for good mapping)
- Fraction of mapped reads
($\geq 90\%$ in good samples)

Paired-end libraries


- Insert size distribution
(distance between F/R read)
- Fraction of reads with a proper pair
($\geq 90\%$ usually)

Reads processing and sample QC


Further clean of sequencing artifacts
and check for sample contamination

- duplicated reads
- UMI decomposition
- sex check
- contamination estimates

PCR and optical duplicates




PCR or optical duplicates



Duplicated reads


A pile of reads that actually represent multiple identical copies originated from the same DNA molecule!

They can be identified as group of reads that have identical start/end points. They are marked and ignored during variant calling. Available tools: [picard](#), [samblaster](#)




Duplicated reads may affect genotyping when the duplicated molecule contains an error / SNVs

UMI based read deduplication to increase accuracy



A scatter plot showing the relationship between 'Scaled mean depth on X' (X-axis) and 'Scaled mean depth on Y' (Y-axis). The X-axis ranges from 1.0 to 2.0, and the Y-axis ranges from 0.0 to 1.2. Data points are colored by sex: blue for 'male' and red for 'female'. Male fish are clustered in the upper left region, with X values between 1.0 and 1.1 and Y values between 0.75 and 1.15. Female fish are clustered in the lower right region, with X values between 1.9 and 2.0 and Y values near 0.0.

Sex	Scaled mean depth on X	Scaled mean depth on Y
male	1.02	0.75
male	1.03	0.95
male	1.04	0.98
male	1.05	1.00
male	1.06	1.02
male	1.07	1.05
male	1.08	1.08
male	1.09	1.10
male	1.10	1.12
male	1.12	1.15
female	1.92	0.00
female	1.95	0.00
female	1.98	0.00
female	2.00	0.00
female	2.02	0.00




autosomes


Male
Haploid X
scaled depth = 1
Haploid Y
scaled depth = 1

Male
Haploid X
scaled depth = 1
Haploid Y
scaled depth = 1


$$\frac{\text{mean depth chrX}}{\text{mean depth autosomes}} \times 2$$



Additional sample-level QC - detect contamination



Increased fraction of “unexpected” alleles in homozygous genotypes



Increased N of heterozygous genotypes
Increased het GT with outlier AR